

ISO 9001:2008 Co.



# Rehsiff Scientific®

Manufacture & Impotor of  
Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

## AMMONIUM BIFLUORIDE PURE

<b>PRODUCT CODE</b>	<b>R01035</b>
<b>SYNONYMS</b>	Ammonium Hydrogen difluoride
<b>C.I. NO.</b>	—
<b>CASR NO.</b>	1341-49-7
<b>ATOMIC OR MOLECULAR FORMULA</b>	NH <sub>4</sub> F.HF
<b>ATOMIC OR MOLECULAR WEIGHT</b>	57.04
<b>PROPERTIES</b>	—
<b>PARAMETER</b>	<b>LIMIT</b>
Description	White to off white crystalline powder.
Solubility	Soluble in water.
Minimum assay (Acidimetric )	98.0%
<b>Note(s) : Assay (if applicable) method mentioned.</b>	
<b>DANGER</b> <b>Hazard statements</b> :Toxic if swallowed. Harmful if inhaled. May be corrosive to metals. Causes severe skin burns and eye damage. <b>Precautionary statements</b> <b>Prevention:</b> Wash hands thoroughly after handling. Do not eat, drink or smoke when using this product.Do not breathe dust or mist. Wash thoroughly after handling. Use only outdoors or in a well ventilated area. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/clothing and eye/face protection. <b>Response:</b> Call a POISON CENTER or doctor/physician if you feel unwell. IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician. If on skin or hair: remove/take off immediately all contaminated clothing. Rinse with water /shower. Absorb spillage to prevent material damage. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Specific treatment: refer to Label or MSDS. Wash contaminated clothing before reuse. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. <b>Disposal:</b> Dissolve the chemical to be disposed, in water and allow it to run to waste, diluting with large quantitiesof water.	IMDG Code : 8/II UN No. 1727 IATA 8

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: [rehsiffscientific@gmail.com](mailto:rehsiffscientific@gmail.com), [www.rehsiffscientific.com](http://www.rehsiffscientific.com)



ISO 9001:2008 Co.



# Rehsiff Scientific®

Manufacture & Impotor of  
Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

Hazard Pictogram(s) :-



GHS05



GHS06

Replace Date 27-Apr.-2023

**: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)**

Email: [rehsiffscientific@gmail.com](mailto:rehsiffscientific@gmail.com), [www.rehsiffscientific.com](http://www.rehsiffscientific.com)

