

ISO 9001:2008 Co.



# Rehsiff Scientific®

Manufacture & Importer of  
Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

## BISPYRAZOLONE AR (FOR DETERMINATION OF CYANIDE)

<b>PRODUCT CODE</b>	<b>R08058</b>	
<b>SYNONYMS</b>	[3.3-Dimethyl-1,1-Diphenyl-(4,4-bi-2-Pyrazoline)-5,5-Dione]	
<b>C.I. NO.</b>	--	
<b>CASR NO.</b>	(7477-67-0)	
<b>ATOMIC OR MOLECULAR FORMULA</b>	C <sub>20</sub> H <sub>18</sub> N <sub>4</sub> O <sub>2</sub>	
<b>ATOMIC OR MOLECULAR WEIGHT</b>	346.38	
<b>PROPERTIES</b>	--	
<b>PARAMETER</b>	<b>LIMIT</b>	
Description	Light yellow powder.	
Solubility	5% solution in formic acid is clear and yellow.	
Minimum assay (TLC)	98.0%	
<b>Note(s): Assay (if applicable) method mentioned.</b>		
<b>WARNING</b>	IMDG Code : --	
<b>Hazard statements :</b> May cause respiratory irritation. Causes skin irritation. Causes serious eye irritation.	UN No. : --	
<b>Precautionary statements</b>	IATA :--	
<b>Prevention:</b> Avoid breathing dust/fume/gas/mist/vapours/spray. Wash thoroughly after handling. Use only outdoors or in a well-ventilated area. Wear protective gloves/protective clothing/eye protection/face protection.		
<b>Response:</b> IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER or doctor/physician if you feel unwell. If eye irritation persists: Get medical advice/attention.		
<b>Disposal:</b> Dissolve the chemical to be disposed, in water and allow it to run to waste, diluting with large quantities of water. The quantities greater than 10g should be dissolved in water and transferred to heavy metal waste drums for collection by specialist disposal company.		
Hazard Pictogram(s) :		

Replace Date 1 April 2023

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: [rehsiffscientific@gmail.com](mailto:rehsiffscientific@gmail.com), [www.rehsiffscientific.com](http://www.rehsiffscientific.com)

