

# **Rehsiff Scientific<sup>®</sup>**

Manufacture & Importer of

Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

CITRIC ACID ANHYDROUS AR			
Meets Analytical Specification of IP, BP, USP, Ph. Eur.			
PRODUCT CODE	R03179		
SYNONYMS	-		
с.і. no.			
CASR NO.	(77-92-9)		
ATOMIC OR MOLECULAR FORMULA	C <sub>6</sub> H <sub>8</sub> O <sub>7</sub>		
ATOMIC OR MOLECULAR WEIGHT	192.13		
PROPERTIES	Strongly acid test; Very soluble in water and alcohol , soluble in ether. Combustible.		
PARAMETER	LIMIT		
Description	A white crystalline powder, Colourless crystals or granules.		
Solubility	10% solution in water is clear and colourless.		
Identification (By FTIR)	Conforms.		
Assay (Acidimetric; calc onanhyd. Sub.)	99.5 - 101.0%		
MAXIMUM LIMIT OF IMPURITIES			
Residue after ignition	0.02%		
Water	0.5%		
Insoluble matter	0.005%		
Chloride (Cl)	0.001%		
Phosphate (PO <sub>4</sub> )	0.001%		
Sulphate (SO <sub>4</sub> )	0.002%		
Arsenic (As)	0.0001%		
Barium (Ba)	Passes test.		
Calcium (Ca)	0.02%		
Iron (Fe)	0.0003%		
Lead (Pb)	0.0002%		
Heavy metals (as Pb)	0.001%		
Oxalate (C <sub>2</sub> O <sub>4</sub> )	0.05%		
Oxalic acid	0.035%		
Organic volatile impurities	Meets the requirement.		
Readily carbonisablesubstance	Passes test.		
Note(s) : Assay (if applicable) method mentioned 2] Store at : <30°C			

: Regd Office, 201, MM Nagar, Ahmedabad (Gujarat)








## **Rehsiff Scientific<sup>®</sup>**

#### Manufacture & Importer of

Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

<ul> <li>Hazard statements: May cause respiratory irritation. May be harmful if swallowed.Causes mild skin irritation.</li> <li>Causes serious eye irritation. Toxic to aquatic life.</li> <li>Precautionary statements</li> <li>Prevention: Wash hands thoroughly after handling.</li> <li>Response: If skin irritation occurs, seek medical advice/attention. Wear eye/faceprotection. If eye irritation</li> </ul>	IMDG Code : UN No. : IATA :
persists, get medical advice/attention. IF SWALLOWED:Call a POISON CENTER or doctor/physician if you feel unwell. IF IN EYES: Rinsecautiously with water for several minutes. Remove contact lenses, if present and easy	







## **Rehsiff Scientific**<sup>®</sup>

### Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

to do. Continue rinsing.		
<b>Disposal</b> :Dispose of contents and containe	er in accordance with relevant legislation.	
Hazard Pictogram(s) :	GHS05 GHS07	

Replace Date 1 April 2023

: Regd Office, 201, MM Nagar, Ahmedabad (Gujarat) Email: <u>rehsiffscientific@gmail.com</u>, <u>www.rehsiffscientific.com</u>

