

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

| PRODUCT CODE | R08188 | | | |
|--|--|--|--------------------------------------|--|
| SYNONYMS | (TriacetoxyPeriodine, 1,1,1-Tris (Acetyloxy)(1,1- Dihydro-1,2-Benziodoxol-3(1H) | | | |
| C.I. NO. | | | | |
| CASR NO. | (87413-09-0) | | o #0 | |
| ATOMIC OR MOLECULAR FORMULA | C ₁₃ H ₁₃ IO ₈ | [[o > ch | | |
| ATOMIC OR MOLECULAR WEIGHT | 424.14 | | 0/6-4 | |
| | | | H ₃ C → O CH ₃ | |
| PROPERTIES | | | 1000 | |
| PARAMETER | LIMIT | | | |
| Description | White powder, Chips, crystalline powder and/or chunks. | White powder, Chips, crystalline powder and/or chunks. | | |
| Solubility | Soluble in acetone. | | | |
| Minimum assay | 97.0% | | | |
| Melting range | 130 – 133°C | | | |
| Note(s): Assay (if applicable) method me | ntioned. | | | |
| DANGER Hazard statements: May cause respiratory irritation. May intensify fire; oxidizer. May be harmful if swallowed. Causes mild skin irritation. Causes serious eye irritation. | | IMDG Co | de : | |
| | | 5.1 / IIIUN | | |
| Precautionary statements | | No. | 1479 | |
| Prevention: Wash hands thoroughly after handling. Take any precaution to avoid mixing with combustible or | | IATA | : 5.1 | |
| ncompatible materials. Keep away from heat. | | | . 3.1 | |
| | advice/attention. Wear eye/face protection. If eye irritation persists, alla POISON CENTER or doctor/physician if you feel unwell. IF IN EYES: | | | |
| | Remove contact lenses, if present and easy to do. | | | |
| Continue rinsing. | remove contact tenses, if present and easy to do. | | | |
| Disposal: The quantities greater than 10g should | be dissolved in dilute acid and transferred to heavy metal waste drum | for collection | on by specialist | |
| lisposal company. | | | | |

replace Date 1 April 2023

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com

