

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

DIETHYL ETHER AR			
PRODUCT CODE	R05900		
SYNONYMS	Solvent ether, Ethyl ether		
C.I. NO.	N/A		
CASR NO.	60-29-7		
ATOMIC OR MOLECULAR FORMULA			
	(C₂H₅)₂O 74.12		
ATOMIC OR MOLECULAR WEIGHT	Caution: Although, this material is not		
PROPERTIES	anhydrous, Diethyl ether tends to form explosive peroxides, especially when anhydrous. It should not be allowed to evaporate to dryness unless the absence of peroxides has been confirmed.	^ ^	
	Hydroscopic with aromatic odor Storage Temp. : 430°C	₃ C´ `O´ `CH ₃	
PARAMETER	LIMIT		
Description	A clear colourless liquid with a characteristic odour having o	olour not more than 10 Hazen units.	
Solubility	Miscible with alcohol.		
Minimum assay(By GC) Refractive index n20	99.5%		
D D	1.351-1.353		
Wt. per ml at 20°C	0.713-0.715 g		
Boiling range (34-35°C)	95%		
MAXIMUM LIMIT OF IMPURITIES			
Acidity(CH₃COOH)	0.001%		
Water	0.20%		
Residue on evaporation (excluding stabiliser) Aldehydes & Ketones	0.001% 0.002%		
Ethanol (C ₂ H ₅ OH)	0.05%		
Methanol (CH ₃ OH)	0.05%		
Peroxide (H ₂ O ₂)	0.00003%		
Substances darkened by Sulphuric acid	Passes test		
Corbonyl compound (CO)	0.001%		
Stabiliser BHT	Abt 10 ppm		
Note(s):Assay (if applicable) method mention	ed.		
Causes skin irritation. Causes eye irritation. May cause PRECAUTIONARY STATEMENTS: Prevention: Do not eat, drink or smoke when using to non-sparking tools. Ground/bond container and rece Contaminated clothing should not be allowed out of a Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves and eye/faceprotection. Keep againststatic discharge. Keep away from heat/sparks, Response: IF SWALLOWED: Rinse mouth. Do NOT inconserveral minutes. Remove contact lenses, if present an before reuse. If skin irritation or rash occurs, seek medical soap and water. Remove/Take off immediately all compared to the serveral minutes.	his product. Wash hands thoroughly after handling. Use on iving equipment. Wash thoroughly after handling. the workplace. Use explosion-proofelectrical/ventilating/lighting/equipment of container tightly closed. Take precautionary measures for one flame - No smoking. duce vomiting. IF IN EYES: Rinse cautiously with water for and easy to do. Continue rinsing. Wash contaminated clothing edical advice/attention. If eye irritation persists, get medical advice/attention. IF ON SKIN: Gently wash with plenty of intaminated clothing. Specific treatment: refer to Label or the reuse. If onskin or hair: remove/take off immediately all	IATA 3 ent.	

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com





Rehsiff Scientific®

Manufacture & Importer of

Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

Disposal: Allow to evaporate in safe open area (applicable for small quantities). The chemical should be burnt in a chemical incinerator equipped with burner and scrubber. If the volumes are larger than 500ml and no incinerator

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com





Rehsiff Scientific®

Manufacture & Importer of

Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

and site of evaporation is not available then transfer the waste chemical to waste solvent drum for collection by specialist disposal company. After testing for peroxide the product should be given to specialized disposal	Plagnostic Nit & Eds We
company. Hazard Pictogram(s):- GHS05 GHS07	

Replace Date 1 April 2023

GHS02

GHS07

 $\textbf{Email:} \ \underline{\textbf{rehsiffscientific@gmail.com}}, \underline{\textbf{www.rehsiffscientific.com}}$

