## ISO 9001:2008 Co.

## Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

DIETHYLENE TRIAMINEFOR SYNTHESIS		
PRODUCT CODE	R09225	
	Bis(2-aminoethyl)amine, 2,2'-	_
SYNONYMS	Diaminodiethylamine  Diaminodiethylamine	
C.I. NO.	N/A LL N	N⊔.
CASR NO.	111-40-0	N N N N N N N N N N N N N N N N N N N
ATOMIC OR MOLECULARFORMULA	C4H13N3	Н
ATOMIC OR MOLECULARWEIGHT	103.17	
PROPERTIES		
PARAMETER	LIMIT	
Description	A yellow viscous liquid, hygroscopic having ammonical odour.	
Solubility	Miscible with water forming clear solution.	
Minimum assay(By GC)	97.0%	
Wt. per ml at 20°C	0.947 - 0.950 g	
Refractive Index(n) 20 D	1.483 - 1.485 g	
MAXIMUM LIMIT OF IMPURITY		
Water	1.0%	
Note(s):Assay (if applicable) method men	ntioned.	
allergic skin reaction. May be corrosive to metalife.  PRECAUTIONARY STATEMENTS:  Prevention: Do not eat, drink or smoke when u dust or mist. Wash hands thoroughly after hand Contaminated clothing should not be allowed o ust/fume/gas/mist/vapours/spray. Wash thorouprotection.  Response: Wash contaminated clothing before and keep at rest in a position comfortable for be CENTER ordoctor/physician. If skin irritation or a cautiously with water for several minutes. Remo SWALLOWED: Rinse mouth. Do NOT induce von position comfortable for breathing. Specific treatments. If on skin or hair: remove/take off immediates.	reuse. IF INHALED: If breathing is difficult, remove to fresh air reathing. If experiencing respiratory symptoms call a POISON rash occurs, seek medical advice/attention. IF IN EYES: Rinse ove contact lenses, if present and easy to do. Continue rinsing. IF niting. IF INHALED: Remove to fresh air and keep atrest in a stament: refer to Label or liately all contaminated clothing. Rinse withwater/shower. Absorb ON CENTER or doctor/physician if you feel unwell.IF ON SKIN:	
	ombustible solvent and burnt in a chemical incinerator quipped	-
with burner and scubber.		

## : Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: <a href="mailto:rehsiffscientific@gmail.com">rehsiffscientific.com</a>, <a href="mailto:www.rehsiffscientific.com">www.rehsiffscientific.com</a>





## Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

Replace Date 1 April 2023

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

 $\textbf{Email:} \ \underline{rehsiffscientific@gmail.com}, \underline{www.rehsiffscientific.com}$ 

