

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

	L (+) TARTARIC ACID AR	
PRODUCT CODE	R02680	
SYNONYMS	(dextro-rotatory (+) tartaric acid)	
C.I. NO.		
CASR NO.	(87-69-4)	
ATOMIC OR MOLECULAR FORMULA	(CHOH.COOH)₂	
ATOMIC OR MOLECULAR WEIGHT	150.09	
PROPERTIES	Odourless, acidic taste and optically active	
PARAMETER	LIMIT	
Description Solubility	Colourless crystals or white or almost white crystalline por 10% solution in water is clear and colourless to not more t	
Identification	Passes test.	
Assay (on dried sub.; Acidimetric)	99.5 – 100.5%	
Specific Rotation [α] $_{D}^{20}$ (C=20%soln in water)	+ 12°to + 12.8°	
MAXIMUM LIMIT OF IMPURITIES		
Loss on drying (at 105°C)	0.2%	
Sulphated ash	0.01%	
Insoluble matter	0.005%	
Chloride (Cl)	0.0005%	
Phosphate (PO ₄)	0.0005%	
Sulphate (SO ₄)	0.01%	
Arsenic (As)	0.00002%	
Copper (Cu)	0.00005%	
Calcium (Ca)	0.002%	
lron (Fe)	0.0003%	
Heavy metals (as Pb)	0.0005%	
Lead (Pb)	0.0005%	
Magnesium (Mg)	0.002%	
Oxalate (C ₂ O ₄)	0.01%	
Organic volatile impurities	Meets the requirements.	
Note(s): 1) Assay (if applicable) method mer 2) "Being natural products pale yellow colo		
WARNING		IMDG Code:
Hazard statements: Causes eye irritation.		UN No. :
Precautionary statements Prevention: Wash hands thoroughly after handling.		IATA :
Response: If eye irritation persists, get medical adv several minutes. Remove contact lenses, if presen	rice/attention. IF IN EYES: Rinse cautiously with water for tand easy to do. Continue rinsing.	

: Regd Office, 201, MM Nagar, Ahmedabad (Gujarat)

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com





Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

Disposal Add in small quantities to large, stirred excess of water, keeping the final concentration less than 2%. Neutralize with 5% sodium hydroxide soln. and run to waste with large volumes of running water. The chemical should be mixed combustible solvent and burnt in a chemical incinerator equipped with burner and scrubber..



Replace Date 1 April 2023

GMP



Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

ŀ	Hazard pictogram (s) :	

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com

