

## Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

PRODUCT CODE	R08297	
SYNONYMS		_
C.I. NO.		
CASR NO.	(7447-41-8)	
ATOMIC OR MOLECULAR FORMULA	LiCl	
ATOMIC OR MOLECULAR WEIGHT	42.39	LiCl
PROPERTIES	Hygroscopic. Not to be used as dietary salt substitute.	
PARAMETER	LIMIT	
Description	A white hygroscopic granular / crystalline powder.	
Solubility	Soluble in water	
Minimum Assay(Argentometric)	99.0%	
pH (of 5% Solution in water)	5-7	
MAXIMUM LIMIT OF IMPURITIES		
Loss on drying	1.0%	
Insoluble matter	0.01%	
Titrable base	0.008 meq /g	
Nitrate ( NO <sub>3</sub> )	0.001%	
Sulphate (SO <sub>4</sub> )	0.005%	
Barium (Ba)	0.002%	
Calcium (Ca)	0.005%	
Heavy metal (as Pb)	0.002%	
Iron (Fe)	0.0005%	
Lead (Pb)	0.0005%	
Magnesium (Mg)	0.005%	
Potassium (K)	0.01%	
Sodium (Na)	0.2%	
Note(s) : Assay (if applicable) method ment		
DANGER		IMDG Code :
	auses skin irritation. Causes eye irritation.May damage the unborn child	UN No. :
Precautionary statements		IATA :
	ecautions have been read andunderstood. Use personal protective equipment Obtain special instructions before use. Do not eat, drink or smoke when using	
as required.Wash thoroughly afternandling. this product.Wash hands thoroughly after h	,	
	water for several minutes. Removecontact lenses, if present and easy to do.	
•	et medical attention advice.If eye irritation persists, get medical	
advice/attention. Ifskin irritation occurs, see	k medical advice /attention. Remove /Take off immediately allcontaminated	
<b>.</b>	MSDS.IF ON SKIN: Gentlywash with plenty of soap and	
water.Wash/Decontaminate removed clothin	g before reuse.	

## : Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com





## Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

: Regd Office, 201,MM Nagar, Ahmedabad (Gujarat)

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com

