

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

ACETIC ACID GLACIAL FOR SYNTHESIS (MEETS ANALYTICAL SPECIFICATION OF IP) R01003 SYNONYMS N/A C.I. NO. N/A CASR NO. 64-19-7 CH₃COOH ATOMIC ORMOLECULAR FORMULA 60.05 ATOMIC ORMOLECULAR WEIGHT PROPERTIES Combustible **PARAMETER** LIMIT Description A clear liquid, not more than 10 Hazen units in colour. Solubility Miscible with water forming a clear &colourless solution. 99.5% Minimum assay (acidimetric) Wt.per ml at 20°C 1.048 - 1.051 g Min. 15.8°C Freezing point MAXIMUM LIMIT OF IMPURITIES Non volatile matter 0.01% Chloride (CI) 0.0005% Formic acid (HCOOH) 0.1% Sulphate (SO_4) 0.0005% Arsenic (As) 0.0001% Iron (Fe) 0.0001% Lead (Pb) 0.0001% Reducing substances Passess test Filtered through 0.2 um membrane filter. Packed under inert gas. Note(s): Assay (if applicable) method mentioned. IMDG Code: 3(8)/II DANGER HAZARD STATEMENTS :Flammable liquid and vapour. Harmful in contact with Harmful if inhaled. May be JN No. 2789 harmful if swallowed. May be corrosive to metals. : 8 Causes severe skin burns and eye damage. Causes serious IATA :3 (8) eye damage. Harmful to aquatic life. PRECAUTIONARY STATEMENTS: Prevention: Keep away from heat/sparks/open flames/hot surfaces. – No moking. Keep container tightly closed. Keep only in original container. Ground/bond container and receiving equipment. Use explosion-proof electrical/ventilating/lighting equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Do not breathe dust/fume/gas/mist/vapours/spray. Avoid breathing

: Regd Office, 201, MM Nagar, Ahmedabad (Gujarat)

contaminated clothing before reuse. Absorb spillage to prevent material damage.

dust/fume/gas/mist/vapours/spray. Wash thoroughly after handling. Use only outdoors or in a well-ventilated area. Avoid release to the environment. Wear protective gloves/protective clothing/eye protection/face protection. Response: IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF ON SKIN: Wash with plenty of soap and water. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. IF INHALED: Remove to fresh air and keep at rest in a position comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician if you feel unwell. Wash

Email: rehsiffscientific@gmail.com, www.rehsiffscientific.com





Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

Disposal: Add in small quantities to large, stirred excess of water, keeping the final concentration less than 2%. Neutralize with 5% sodium hydroxide soln. and run to waste with large quantities of running water. Hazard pictogram(s):		

Replace Date 1 April 2023

