

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

	Laboratory Fine Chem	icals, Diagnostic Kit & Lab w	
ADR	ENALINE HYDROGEN TARTRATE AR		
	Meets IP,BP,USP,Ph,Eur		
PRODUCT CODE	R04040		
SYNONYMS	(Adrenaline hydrogen tartrate, L-Adrenaline-D-	au.	
C.I. NO.	hydrogen tartrate)	HN CH ₃	
ASR NO.	(51.42.2)	НО	
	(51-42-3) C₀H₁₃NO₃.C₄H₅O₅	о он	
ATOMIC OR MOLECULAR FORMULA	333.3	ОН	
TOMIC OR MOLECULAR WEIGHT PROPERTIES	It darkens on exposure to air and light.		
PARAMETER	LIMIT		
Description	White or gravish white crystalline powder. It darkens on e	White or grayish white crystalline powder. It darkens on exposure to air and light.	
olubility	5% solution in water is clear		
dentification	Passes test.		
Minimum assay(Non-agueous)			
oH 1% Aq. solution	99.0% 2.8 – 4.0		
Melting range	2.6 – 4.0 147° – 152°C		
MAXIMUM LIMIT OF IMPURITIES			
oss on drying (on P ₂ O ₅)	0.5%		
Sulphated ash	0.1%		
Phenones /Adrenalone (0.2%in 0.1N HCl at 3	10 nm) 0.1 AU		
Noradrenaline (TLC)	Passes test.		
Note(s): Assay (if applicable) method ment	ioned.		
DANGER		IMDG Code : 6.1/I	
Hazard statements: May cause respiratory irritation. Toxic if inhaled. Toxic in contact with skin. Toxic if swallowed. Causes mild skin irritation.		UN No. 2811	
dust/fume/gas/mist/vapours/spray. Avoid Ismoke when using this product. Wash hands	ng. Use only outdoors or in a well ventilated area. Do not bre breathing dust/fume/gas /mist/vapours/spray. Do not eat, drin thoroughly after handling. ENTER or doctor/physician. If skin irritation occurs, seek me	k or	
advice/attention. Keep container tightly clo comfortable for breathing. IF ON SKIN: Ge	ised. IF INHALED: Remove to fresh air and keep at rest in a posintly wash with plenty of soap and water. Call a POISON CENTE treatment: refer to Label or MSDS. IF SWALLOWED: Immediately of soap and water.	ition R or	
	re/Take off immediately all contaminated clothing.		
	ed, in water and allow it to run to waste, diluting with large quanti	ties	
Hazard Pictogram(s) :-			

Replace Date 1 April 2023

