

## Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

| AJMALINE PURE  |   |   |
|--|---|---|
| PRODUCT CODE   | R08727  | <u></u>   |
| SYNONYMS   | <del></del>   | .OH   |
| C.I. NO.   |   | \(\sigma_1\)                                    |
| CASR NO.   | (4360-12-7)   |   |
| ATOMIC OR MOLECULAR FORMULA  | $C_{20H_{26}N_2O_2}$  | N OH  |
| ATOMIC OR MOLECULAR WEIGHT   | 326.44  | CH <sub>2</sub>                                 |
| PROPERTIES   | <del></del>   | CH <sub>3</sub> CH <sub>2</sub> CH <sub>3</sub> |
|  |   |   |
|  |   |   |
| PARAMETER  | LIMIT   |   |
| Description  | White to pale yellow crystalline powder.                              |   |
| Solubility   | Soluble in water.   |   |
| Minimum assay  | 98.0%   |   |
| Melting point  | About 189°C   |   |
| Note(s) : Assay (if applicable) method mention                       | oned.   |   |
| WARNING  |   | IMDG Code :                                     |
| , , ,  | ritation. Causes skin irritation. Causes serious eye irritation.      | UN No. :  |
| Precautionary statements Prevention: Avoid breathing dust/fume/gas/i | mist/vapours/spray. Wash thoroughly after handling. Use only          | IATA :  |
|  | rotective gloves/protective clothing/eye protection/face protection   | n.  |
| •  | nd keep at rest in a position comfortable for breathing. IF IN EYES:  |   |
| •  | es. Remove contact lenses, if present and easy to do. Continue        |   |
|  | ician if you feel unwell. If eye irritation persists: Get medical     |   |
| advice/attention.  | d, in water and allow it to run to waste, diluting with large quantit | ios   |
|  | uld be dissolved in water and transferred to heavy metal waste dr     |   |
| for collection by specialist disposal company.                       |   |   |
| Hazard Pictogram(s) :  |   |   |
|  | ^ ^   |   |
|  |   |   |
|  |   |   |

Replace Date 1 April 2023

