

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

| PRODUCT CODE | R08132 | |
|---|---|--|
| SYNONYMS | α-Chloroacetamide | _ |
| C.I. NO. | | 0 |
| CASR NO. | (79-07-2) | . [] |
| ATOMIC OR MOLECULAR FORMULA | C₂H₄CINO | NH _o |
| ATOMIC OR MOLECULAR WEIGHT | 93.51 | 11112 |
| PROPERTIES | Soluble in water and alcohol, insoluble in ether. | |
| PARAMETER | LIMIT | |
| Description | Colourless to pale yellow crystal with characteristic odour. | |
| Solubility | 5% solution in water and alcohol should be clear. | |
| Minimum assay (GC) | 97.0% | |
| Melting range | 115 - 118 ⁰ C | |
| Note(s) : Assay (if applicable) method men | tioned. | |
| prolonged or repeated exposure. Harmful to Precautionary statements Prevention: Do not breathe dust/fum dust/fume/gas/mist/vapours/spray. Use onloroduct. Washhands thoroughly after hand Response: Get medical advice/attention if y Gently wash with plenty of soap and water preathing. Keep container tightly closed. Co | ne/gas/mist/vapours/spray. Wear protectivegloves/clothing Avoid by outdoorsor in a well ventilated area. Do not eat, drink or smoke when uling. Tou feel unwell.Immediately call a POISONCENTER or doctor/physician. IF or r.IFINHALED: Remove to fresh air and keep at rest in a position comfor call a POISON CENTER or doctor/physician if you feelunwell.IF SWAL physician.Specific treatment: refer to Label or MSDS. Wash/ Decontamina | IATA : 6.1 preathing using this ON SKIN: table for LOWED: |
| removed clothingbefore reuse. Remove/Take Disposal: Dispose of contents and container | | |

Replace Date 1 April 2023

