

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

| | DEAE CELLULOSE 52 AR | | |
|--|--|--------|------|
| | FOR COLUMN CHROMATOGRAPHY | | |
| PRODUCT CODE | R09716 | | |
| SYNONYMS | (Diethylaminoethyl-cellulose) | | |
| C.I. NO. | | | |
| CASR NO. | (9013-34-7) | | |
| ATOMIC OR MOLECULAR FORMULA ATOMIC OR MOLECULAR WEIGHT PROPERTIES | | | |
| PARAMETER | LIMIT | | |
| Description | White to off-white powder or solid | | |
| Solubility | Insoluble in water | | |
| Capacity | 0.80 - 1.10 meq/g, dry weight | | |
| MAXIMUM LIMIT OF IMPURITIES | | | |
| ONases | None detected | | |
| RNases | None detected | | |
| Water (KF) | 65.0% | | |
| Note(s): Assay (if applicable) method | mentioned. | | |
| extinction. If on skin or hair: remove/take of HAZARD STATEMENTS: Highly flammable li PRECAUTIONARY STATEMENTS: Prevention: Take precautionary measures a electrical/ventilating/lighting/equipment. G | ENTER or doctor/physician if you feel unwell. In case of fire, use foam for if immediately all contaminated clothing. Rinse with water/shower. quid and vapour. May be harmful if swallowed. against static discharge. Use explosion-proof round/bond container and receiving equipment. Use only non-sparking tools. Io smoking. Keep container tightly closed. Wear protective gloves and | IATA 3 | .170 |
| Hazard Pictogram(s) : | GHS02 | | |

Replace Date 1 April 2023

