

Rehsiff Scientific®

Manufacture & Importer of Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

| 4- METHOXY BENZYLAMINE FOR SYNTHESIS | | |
|---|-----------------------------------|------------------|
| PRODUCT CODE | R09497 | |
| SYNONYMS | | |
| C.I. NO. | | NH ₂ |
| CASR NO. | (2393-23-9) | <u> </u> |
| ATOMIC OR MOLECULAR FORMULA | C ₈ H ₁₁ NO | |
| ATOMIC OR MOLECULAR WEIGHT | 137.18 | |
| PROPERTIES | | Y |
| | | OCH ₃ |

| PARAMETER | LIMIT | |
|------------------------------------|--|--|
| Description | Colourless to yellow liquid material darkens in storage with no loss of purity | |
| Solubility | Soluble in water. | |
| Minimum assay (GC) | 98.0% | |
| Wt. per ml at 25°C | About 1.05 g | |
| Refractive index (n) ²⁰ | About 1.546 | |
| Boiling range | 236-237°C | |

Note(s): Assay (if applicable) method mentioned.

WARNNING

HAZARD DANGER Harmful if swallowed. Causes skin irritation. Causes serious eye damage. May cause respiratory UN No. 2735 irritation.

IATA 8

PRECAUTIONARY STATEMENTS

Prevention: Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray. Wash. thoroughly after handling. Do not eat, drink or smoke when using this product. Use only outdoors or in a well- ventilated area. Wear protective gloves/protective clothing/eye protection/face protection.

Response: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Wash with plenty of soap and water. IF INHALED: Remove victim to fresh air and keep at rest in a position. Comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove. Contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. Call a POISON CENTER or doctor/physician if you feel unwell. Rinse mouth. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash before re- use

Disposal: The chemical should be mixed with combustible solvent and burnt in a chemical incinerator equipped with burner and scrubber

Hazard Pictogram(s):





Replace Date 1 April 2023

