

Rehsiff Scientific[®]

Manufacture & Importer of

Laboratory Fine Chemicals, Diagnostic Kit & Lab wares

| PRODUCT CODE | R08396 | |
|--|---|-----------|
| SYNONYMS | (α-Nitroso-β-naphthol) | NO |
| C.I. NO. | - | NO |
| CASR NO. | (131-91-9) | OH |
| ATOMIC OR MOLECULAR FORMULA | C ₁₀ H ₆ (OH).NO | |
| ATOMIC OR MOLECULAR WEIGHT | 173.17 | |
| PROPERTIES | Soluble in benzene, ether, caustic alkali solutions glacial acetic acid. | |
| PARAMETER | LIMIT | |
| Description | Brown or yellowish brown crystalline powder. | |
| Solubility | Soluble in alcohol, 2% solution in 50% hot acetic acid | is clear. |
| Vinimum assay (dried) | 98.0% | |
| Melting point | 104-108°C(dec) | |
| MAXIMUM LIMIT OF IMPURITIES | | |
| Sulphated ash | 0.2% | |
| Loss on drying (70 ⁰ C) | 3.0% | |
| Sensitivity for CO | To pass the test. | |
| Note(s):Assay (if applicable) method me | entioned. | |
| Precautionary statements Prevention: Avoid breathing dust/fume/gas putdoors or in a well-ventilated area. Wear p Response: IF INHALED: Remove to fresh air a Rinse cautiously with water for several mini- rinsing. Call a POISON CENTER or doctor/physi advice/attention. | tation. Causes skin irritation. Causes serious eye irritation. s/mist/vapours/spray. Wash thoroughly after handling. Use only rotective gloves/protective clothing/eye protection/face protection. and keep at rest in a position comfortable for breathing. IF IN EYES: utes. Remove contact lenses, if present and easy to do. Continue cian if you feel unwell. If eye irritation persists: Get medical | |
| of water. Hazard Pictogram(s) : | | |

Replace Date 1 April 2023

